

## → Who is MyWay suitable for?

MyWay is an upright movement device which supports children who cannot stand or walk independently in an upright position and allows them to move about their environment.

MyWay is intended for use by children already using a supportive walking device; children with Cerebral Palsy in GMFCS levels III–IV; children with Cerebral Palsy in GMFCS level V who have sufficient range of movement in their lower limbs to enable passive/active stepping; and children with a physical disability which prevents them from independent upright mobility.



## → How many sizes of MyWay are there?

MyWay is available in 3 frame sizes with 6 corresponding harnesses.

## → What is the maximum user weight/size?

	Size 1 Frame	Size 2 Frame	Size 3 Frame
Age (approx.)	1-5	4-11	9-16
User Shoulder Height (to floor) (mm)	635 – 920	870 – 1190	1070 – 1490
User Saddle Height (inside leg) (mm)	220 – 440	400 – 610	550 – 810
Max. User Weight (kg)	25	50	80

	Harness S1	Harness S2	Harness S3	Harness S4	Harness S5	Harness S6
Age (approx.)	1 – 4	2 – 6	5 – 9	8 – 12	12 – 15	16 – 18
Waist at navel level (cm)	46 – 56	50 – 60	55 – 65	60 – 70	70 – 80	80 – 90
Max. User Weight (kg)	50	50	50	80	80	80
Colour	Purple	Green	Blue	Red	Orange	Grey

**Please note:** The maximum user weight for a MyWay is determined by the combination of frame and harness sizes selected- the maximum user weight is the lower of these limits. For example, if a Size 3 frame and a Size 3 harness are selected, the maximum user weight is 50kg, as determined by the harness limit.

## → Are there any accessories available?

Yes! There are a range of optional accessories available including:

- **Adult Connection Belt** – this allows an adult to attach to the MyWay frame and control the frame's movements, providing additional support to the child if required. Using the connection belt allows the adult to facilitate walking and therapy exercises whilst leaving their hands free to interact with the child. Please note, the adult handles must also be purchased in order to use the Adult Connection Belt.
- **Saddle**
- **User Handles**
- **Ankle Guides**
- **Headrest**
- **Adult Handles**
- **Hoist Straps**

## → Is MyWay suitable for use indoors and outdoors?

MyWay is suitable for use indoors and outdoors on smooth, level surfaces. MyWay should not be used on uneven or rough terrain. Potential environments for use should always be risk assessed to ensure the child's safety.

## → Does MyWay fold?

Yes! The MyWay Size 2 and Size 3 can be split into two parts and the wheeled frame is foldable. The Size 1 MyWay does not split due to its smaller footprint but it is also foldable.

## → How heavy is MyWay?

	Size 1 Frame	Size 2 Frame	Size 3 Frame
Frame Weight (kg)	12 (S1 frame and upper support unit do not separate)	6.5	9.5
Upper Support Unit Weight (kg)		10.5	11

**→ What are MyWay's dimensions?**

	Size 1 Frame	Size 2 Frame	Size 3 Frame
Frame Width (mm)	570	670	720
Frame Length (mm)	720	850	1100

**→ Is the harness washable?**

The MyWay harness can be machine washed at 40°C on a gentle program. A washbag is provided with your MyWay harness for convenience. Furthermore, most stains can be removed using a clean microfibre cloth and lukewarm water.

**→ Can I use any existing hoist straps I have with the MyWay harness?**

No, only the MyWay hoist straps should be used with the MyWay harness.

**ALLIED  
MEDICAL**  
LIMITED