



**ALLIED
MEDICAL**
LIMITED

USER MANUAL

Swivel Cushion

MP234745



Soft Swivel Transfer Seat



**Please inspect all parts carefully
for cracks/breaks prior to use.**


USER INSTRUCTIONS

This cushion is designed to aid in transfer in and out of automobile or chair.

- Soft top-surface fabric
- Anti-slip bottom
- 16" diameter
- Lightweight for easy transport
- Does not contain latex

CARE INSTRUCTIONS

- Do not immerse in water.
- Do not dry clean.
- Do not iron.
- Wipe clean with damp cloth.



Maximum
User Wt:

