

# Otto Bock®

QUALITY FOR LIFE



## Nurmi Neo – Posterior Walker

- Promotes an active, upright posture and natural movements
- Individually adaptable to different body heights and abilities
- Foldable and easy to transport
- "Grows with the child"

## Nurmi Neo – Natural Walking

Nurmi Neo is a walking aid that is positioned behind the body, opening out into the walking direction. This promotes an upright posture and natural walking movements. Existing residual functions can be used and promoted efficiently. Together with its light weight and extra smooth rolling characteristics of the wheels, Nurmi Neo requires little effort for autonomous movement.

The Nurmi Neo walking aid offers a high degree of security thanks to its advanced stability, its reverse-roll locks and its optional anti-tipper device. Foldable and with a light weight of 5 kg, Nurmi Neo is easy to transport and store.

Nurmi Neo "grows with the child". Height-adjustable grip bars allow it to be easily adjusted to different body heights.

A large range of accessories for different requirements and medical conditions allow for versatile and long-term use.

Nurmi Neo: for natural movements and real mobility.



# Nurmi Neo

## Nurmi Neo in Detail

- 3 sizes
- Swiveling front wheels
- Height and depth adjustable grip bars
- Reverse-roll locks at the rear wheels (can be activated and deactivated)
- Foldable
- Driving brake and wheel lock
- Friction brake
- Caster swivel lock
- Anti-tipper
- Standard grip bars, grip bars with universal grips, grip bars with forearm supports and vertical hand grips or "mini" grip bars
- Fold-up seat
- Basket

## Technical Data

Dimensions (cm) and weight (kg)	Size 1	Size 2	Size 3
Push handle height: Grip bar, standard	45 - 57	58 - 72	64 - 78
Distance betw. grips: Grip bar, standard	41	46	51
Push handle height: Grip bar w. univ. grips	38 - 57	51 - 72	57 - 78
Distance betw. grips: Grip bar w. univ. grips	27 - 41	32 - 46	37 - 51
Push handle height: "Mini" grip bar	31 - 48		
Distance betw. grips: "Mini" grip bar	20 - 41		
Grip diameter	3	3	3
Distance forearm supports – floor	53 - 64	66 - 79	71 - 84
Distance betw. forearm supports (center)	25 - 41	30 - 46	35 - 51
Forearm support width x length (inside)	7 x 16	7 x 21	7 x 21
Overall width	63	68	72
Overall length (with anti-tipper +10 cm)	65	76	83
Folding dimensions (L x W x H)	65 x 63 x 24	76 x 68 x 24	85 x 72 x 24
Front wheel diameter	15	15	15
Rear wheel diameter	15	20	20
Load capacity	25	40	55
Weight	5.5	6.0	6.5
Seat height	30	45	48
Turning diameter (without accessories)	86	99	106



Combined gluteal and pelvic supports



Driving brake and wheel lock



Grip bar with forearm supports and vertical hand grips



Grip bar with universal grips for grip width reduction



Storage basket



Anti-tipper



Reverse-roll lock at the rear wheels (can be activated and deactivated)



Friction brake at the rear wheels

Ottobock®

QUALITY FOR LIFE

29 Triton Drive  
PO Box 302 250  
North Harbour  
Auckland, New Zealand

0800 31 61 81  
Ph: +64 9 415 1685  
Fax: +64 9 415 1686  
www.alliedmedical.co.nz



ALLIED  
MEDICAL  
LIMITED